

Validation of an Artificial Intelligence Left Ventricular Ejection Fraction Quantification Software with Cardiac Magnetic Resonance Imaging in Consecutive Non-Selected Patients




Ashton C. Lai, Solomon Bienstock, Rajeev Samtani, Frans Beerkens, Lori Croft, Steve Liao, Eric Stern, Martin E. Goldman,
Icahn School of Medicine at Mount Sinai, New York, NY

Background

- Current quantification of LVEF is subjective, time-consuming, and error-prone.
- Automated artificial intelligence (AI) programs have been developed to reduce variability and increase throughput.
- LVivoEF by DIA®, new artificial intelligence (AI) software, provides an almost instantaneous LVEF quantification from a single-plane apical 4-chamber view.
- We sought to assess the ability of LVivoEF to track and quantify routine clinical TTEs.

Methods

- Retrospective single center study of 111 non-selected random patients (pts) who underwent both TTE and c-MR within 6 months (mean 28 ± 40.3 days) with no interval significant clinical events
- AI-derived LVEF from a single plane 4-chamber view was compared to c-MR in 111 pts. Using linear regression and chi square, the AI 4-chamber LVEF were compared to the cMRI-derived LVEF.
- High-quality images were defined based on endocardial definition, lack of apical foreshortening, lack of wall "dropout," overall signal-to-noise, and vertical axis by two independent physicians.



Results

- Of the 111 pts (53% men; mean age 54.5 years), 59 had a normal EF (>51%), 16 mildly abnormal (41-51%), 15 moderately abnormal (30-40%), and 21 severely abnormal (<30%) while 29 had wall motion abnormalities.
- The correlation between AI and c-MR was highly significant with $R^2 = 0.759$. When only high-quality images were utilized, the correlation between AI and c-MR was $R^2 = 0.918$.
- Finally, AI was highly sensitive (91.67%) and specific (81.33%) for detecting reduced LVEF $\leq 40\%$ versus LVEF > 40% ($\chi^2 = 53.09$, $p < 0.00001$).

Conclusions

- Using a single plane apical 4-chamber view, LVivoEF by DIA® can quantify LVEF and categorize reduced LVEF regardless of image quality in non-selected pts.
- Thus, AI has the ability to dramatically improve and expedite LVEF quantification in routine clinical practice.